


mediCAD[®] 3D Printing

Planning of a Trochleoplasty via 3D Print




mediCAD[®]
3D PRINTING





Technological progress in visualisation techniques like 3D print support clinics in developing customised surgical solutions for patients. Due to the increasing development of technologies and, thus, more complex surgical interventions, better surgery preparation, digital planning and further qualification of doctors are needed. Especially for pre-operation planning 3D models can help create realistic and customised pathologies of the anatomic regions affected and support at therapy decision, or they help visualise the actual situation to doctors-to-be and give them the possibility to simulate the operation in advance.

3D print of the distal femur to prepare a trochleoplasty

Young female patients suffer quite often from pain in their patella caused by a valgus in combination with an augmented antetorsion of the femur. This was also the case of a young patient who went into this clinic. We ask the attending doctor about his experience.



mediCAD® 3D Knee



mediCAD®
3D PRINTING



In this case, the patient has a healthy trochlea at her right leg and a valgus malposition in combination with antetorsion of the femur at her left leg. Additionally, the patient was treated with two failed trochleoplasties which caused a necrosis in this area. Therefore, the attending doctor had the idea to take the healthy side as a model for the pathological side in order to get an adequate possibility to prepare the trochleoplasty. He ordered a 3D print for both trochleae to get a comparison of both and to simulate the operation at the affected trochlea.

Case Study

mediCAD®: Which treatment options do you see for this young female patient?

Doctor: The patient is very young and we have to very carefully reflect and plan what we can “do to her” surgically in this young age.

mediCAD®: Previous conservative therapies have failed?

Doctor: Yes, and there were two failed trochleoplasties because osseous deformities like a pathological TT-TG, a valgus, especially at the left leg, haven’t been considered. Instead of straightening the leg axis, the therapy was reduced to the trochlea hole.

mediCAD®: What can be done better in the case of such a young patient?

Doctor: Generally, for this kind of pain at the patella the range of treatments includes MPFL-reconstruction, tuberosity transfer, trochleoplasty and osteotomies to correct the leg axis and the torsion. Often, we choose combined interventions. In this case, we decided on a medial opening wedge

HTO at the left side to get symmetry of both legs, thus changing also the patella tracking by cutting above the starting point of the patella sinew, the tuberosity tibiae. We start always with interventions of axis correction before pursuing changes of the trochlea or the ligaments.

The screenshot shows the '3D Druck Bestellung' (3D Print Order) interface. At the top, there are steps: '1 STL hochladen' (upload STL), '2 Preis' (price), and '3 Zusammenfassung' (summary). A tip box asks 'Wussten Sie schon?' (Did you know?) and suggests dragging and dropping STL files. Below is an 'UPLOAD' button with a 3D printer icon. The main area is a table with columns: 'Produkt', 'Härtegrad', 'Preis', 'Stückzahl', and 'Summe'. Two items are listed:

Produkt	Härtegrad	Preis	Stückzahl	Summe
 Trochlea_rechts_gesund.stl 3D Druck	Hard	€375.92	1	€375.92
 Trochlea_links_Frakturseite.stl 3D Druck	Hard	€363.06	1	€363.06
Summe				€738.98

services.mediCAD.cloud

For this purpose, he used the mediCAD® 3D Knee Software to segment the area of a CT which he wanted to 3D print and exported the model as .stl file.

Directly from the mediCAD® user interface he opened the mediCAD® service portal for the order process. At this portal he received a price quote for the 3D print. Then everything went really quickly: The doctor entered his delivery address and the print was at his desk only 4 work days later.

mediCAD®: That sounds fascinating. Was this enough to release the patient free of pain from the hospital?

Doctor: We had corrected only one leg and wanted to see first how she was doing with it - even though we knew already that the trochleae was very dysplastic, or was affected by two trochleoplasties which were not working satisfactorily. The decision for the intervention had been taken at a follow-up examination. Most importantly, we now had to do everything right and plan the intervention in detail. Therefore, I wanted to use two 3D prints because I believe that trochleoplasty can be planned very well with this tool.

mediCAD®: Perfect, the future will certainly offer new opportunities for the planning of a trochleoplasty.

Doctor: I agree! In my eyes, 3D planning in combination with the possibility of 3D print will be the future of surgery planning.

mediCAD®: Thank you very much for the interview and your time!



Production Flow

Order Process

1. Register or log in at services.mediCAD.cloud in order to start the order process
2. Upload your saved .stl file of the anatomic region
3. Determine the desired hardness
4. Costs and supply

Production

In collaboration with a mediCAD® 3D Printing partner your order will be processed immediately after order entry. Production is carried out via PolyJet 3D Print.

Delivery


We guarantee an express delivery of your product within 2-4 work days.

Advantages

- Realistic illustration of pathologies
- Communication with the patient
- Surgery preparation and planning
- Anatomic anomalies (deformities, trauma etc.) can be very well illustrated.
- Each bone or affected area can be printed. All you need is a 3D record.
- There are two different degrees of hardness. Depending on the intended use you can select the adequate material, e.g. soft for the reposition at traumas or hard for a realistic simulation of bone hardness.



mediCAD®
3D PRINTING





Digitalisation in clinics and surgeries

Future via new technologies:
3D Print, Web/Cloud, Virtual Reality and Artificial Intelligence

mediCAD Hectec GmbH
Opalstr. 54
DE-84032 Altdorf
Tel. +49 871 330 203-0
E-Mail: info@mediCAD.eu

mediCAD[®]
The Orthopedic Solution
www.mediCAD.eu